



Tall Locking 6 Drawer Cart

Part No. 51033

TrippNT's TL 6D Cart stores small or bulky items. The unique tall slender design creates maximum storage in a small space. This cart features 6 full-extension drawers that all lock with the turn of a key. Heavy-duty, one piece construction adds strength and security for trauma, ER, surgical, anesthesia, critical care, and any hospital, lab or care center use. Perfect for labs, clinics, operating rooms, ER's, patient supplies, hospitals, doctor's offices and similar industries to add extra space, storage, and mobility. Available in 6 color variants. Ships Fully Assembled.



51033



51355



51089





51097



51382




51096




Small Drawer ID Dimensions

W - 11 1/2"
D - 15 3/8"
H - 3 3/8"



Large ID Dimensions

W - 11 1/2"
D - 15 3/8"
H - 7 1/2"



Features:

- 2 Locking & 2 Swivel 4in Twin Wheel Casters, 6 Total Locking Drawers, Top Tray
- Unique Slender Design
- Corrosion-Resistant, HDPE Construction
- Organizes, Stores, Secures, and Mobilizes Supplies



Ships Next Day



Made in the U.S.A



Corrosion Resistant



100% Lifetime Guarantee

Rev: 51033.25