




Core Pipette Fume Hood Cart

Part No. 53602

TrippNT's Core Pipette and Fume Hood Supply Cart with white drawers and pipette slots organizes, stores, and mobilizes pipettes and fume hood supplies. It features a low profile design for storing under fume hoods and biosafety cabinets.



Features:

- 2 Locking & 2 Swivel 3in Casters, 2 Full-Extension Drawers, 6 Pipette Slots
- Designed to store under fume hood or biosafety cabinet
- Corrosion-Resistant, Non-Metal HDPE Construction
- Organizes, Stores, and Mobilizes Pipette Supplies



Ships Next Day



Made in the U.S.A



Corrosion Resistant



100% Lifetime Guarantee

Rev: 53602.1