



24 inch Adjustable Surgical Helmet Rack



Part No. 53941

TrippNT's 24 inch Adjustable Surgical Helmet Rack stores and protects surgical helmets from damage and contamination from falling. Designed in conjunction with a surgical department, this storage rack boasts 3 totally adjustable hooks giving you freedom to optimize placement of your helmets (Styrker® or similar). Also, a great option for hanging other oddly shaped equipment and supplies. A larger 48 inch (5 hook) size is available as well as individual replacement hooks.



Features:

- 24 inch wall mounted base with 3 hooks
- Easy glide adjustable hooks that hand tightens with thumb screws
- Made from white HDPE
- Organizes and protects surgical helmets
- 50 pound weight capacity per hook



Rev: 53941.1



Ships Next Day



Made in the U.S.A



Corrosion Resistant



100% Lifetime Guarantee