




# Height Adjustable Suture Stand Cart

Part No. 53944

The Height Adjustable Suture Stand Cart holds up to 48 suture boxes, with double-sided storage and six slanted shelves. Adjustable from 29 1/4 to 42 1/2" with easy access and smooth mobility. Has two ergonomic built in handles for effortless mobility.



Rev: 53944.1

## Features:

- Provides reliable storage for operating rooms, urgent care, and other medical settings.
- Double-sided with 6 slanted shelves, holds 48 standard suture boxes.
- Made of durable and easy to clean polyethelene and PVC
- Footprint: 24 x 61 1/2 x 19 1/2 inches (WHD)



**Ships Next Day**



**Made in the U.S.A**



**Corrosion Resistant**



**100% Lifetime Guarentee**