



BD Sharps Container 5.4 Quart

Part No. 53293





Protect people and maintain compliance with this durable, disposable sharps container, expertly designed to securely hold needles, syringes, lancets, razors, and other sharp objects. Constructed from high-impact, shatter-resistant polypropylene, its puncture-proof build ensures reliability in demanding environments. The innovative automatic rotating chamber allows for touchless sharps disposal, reducing accidental contact risks and ensuring hygienic safety. Perfect for warehouses, medical facilities, schools, and rest areas, this lightweight container measures 12 1/8" x 10 7/8" x 4 3/4" (WHD) and weighs just 1 pound.

Rev: 53293.1



Features:

- *Touchless disposal with automatic rotating chamber*
- *Puncture-Resistant Design*
- *Securely contains needles, syringes, lancets, razors, and other sharps.*
- *Fits TrippNT Phelbotomy Pro (51536) and Phelbotomy Pro Mini (53292)*



Product Overview



Short Lead Times



Made in the U.S.A



Corrosion Resistant



100% Lifetime Guarantee