



**HEMOBILIZE** © By *Tripp*NT

Softgoods at Bedside


HEMOBILIZE®

## SG Transport Cart™

- Holds up to 60 Softgoods Setups (3 softgoods setups/tray)
- Replenishment system for SG Access Cart™, SG Supply Cart™, and SG Donor Bed™
- Flip lock to easily secure trays
- Incorporated Footbrake
- **Lifetime Warranty on Trays**  
2 Year Warranty on Aluminum Frame



**91103**



### Additional Info

Materials:  
Powder Coated Aluminum  
20 Polycarbonate Trays

Weight: 72 lbs

Dimensions:  
28 1/2" x 68" x 28 1/2"  
Width Height Depth

for more information, visit  
[www.HEMOBILIZE.com](http://www.HEMOBILIZE.com)

Sales & Product Information  
Bob Wright, Plasma Product Specialist  
248.875.1954 | [bobwright@trippnt.com](mailto:bobwright@trippnt.com)



MADE IN USA